


The art of innovation



#### M-24T-Pellet & Film Press.

This press has unique features when it comes to the production of high quality Pellets and Films



#### M-24T-Pellet & Film Press.

All purpose Pellet and Film press for heavy duty applications

Best bench top Press for the production of all sorts of pellets and films. The system can be provided with various accessories dedicated to the application of the client.



#### M-24T-Pellet & Film Press.

M24T-P&F desktop dry manual pellet presses can press up to 24 metric tons. They are frequently used in materials and chemistry research laboratories to prepare compound samples. These presses have built in hydraulic pump with auto return function and 10mm thick transparent safety shield.

#### Specifications


Maximum pressure	24 metric tons (52,910 Lb)
Maximum working space	4x4x5.5" (LxWxH)
Largest compatible die set	2" ID round die sets (2x2" pressing dies will NOT fit)
Hydraulic cylinder travel distance	1 inch
Machine dimension	10.5x8x19" (LxWxH)
Machine weight	84 Lb
Shipping weight	110 Lb
Shipping dimensions (LxWxH)	15x15x22"
Pressing dies	Various diameter dies available at extra cost

Standard package Part description	Quantity (PC)	Part Image
M-24T-Pellet & Film Press.	1	
Oil Tin	1	
Hex L-Key	2	
Double open-end wrench	3	
Sealing o-ring	4	
Operation manual	1	

#### More information!

Ask our sales representatives for more information.


[info@microptik.eu](mailto:info@microptik.eu)



The art of innovation



**13-25mm ID 250°C Heated Dry Pressing Die Set 220V 50/60Hz for Pellet and film production.**  
A set of 10 till - 25mm diameter (ID) hardened steel dry pressing die set for pellet presses, ready for you to make from 10 mm to 25mm diameter cylinder-shape pellets. Heating system with programmable digital temperature controller allows heating up to 250°C.



## M-24T-Pellet & Film Press. 13mm ID 250°C Heated Dry Pressing Die Set 220V 50/60Hz

The M-24T-Pellet & Film Press. Can be provided with various Die sets, depending on the application. We provide also Heated dry pressing die sets for pellets and films. The system comes with a die set et and a temperature controller.



**M-24T-Pellet & Film Press (heated Dry Pressing Set with temperature controller) .**

The controller provides accurate temperature control for the heating of the dies. If required a full software control option can be chosen.



### Specifications

Heated Dry Pressing set	Several diameters possible from 10 till 25mm
Temperature range	From ambient to 250°C
Ejecting (pellet / film)accessories	To facilitate the exit of pellets of films after production
Additional aluminum and Teflon sheets available	To facilitate release of pellets and film

Standard package Part description	Quantity (PC)	Part Image
13mm ID 250°C Heated Dry Pressing Die Set 220V 50/60Hz	1	
Controller	1	
Operation manual	2	

### More information!

Ask our sales representatives for more information.

[info@microptik.eu](mailto:info@microptik.eu)